



YPEL5 blocking peptide (CDBP6448)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	YPEL5 yippee-like 5 (Drosophila) [Homo sapiens (human)]
Official Symbol	YPEL5
Synonyms	YPEL5; yippee-like 5 (Drosophila); CGI-127; protein yippee-like 5
Entrez Gene ID	51646
mRNA Refseq	NM_001127399
Protein Refseq	NP_001120871
UniProt ID	P62699