



ZBED2 blocking peptide (CDBP6450)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZBED2 zinc finger, BED-type containing 2 [Homo sapiens (human)]
Official Symbol	ZBED2
Synonyms	ZBED2; zinc finger, BED-type containing 2; zinc finger BED domain-containing protein 2; zinc finger, BED domain containing 2
Entrez Gene ID	79413
mRNA Refseq	NM_024508
Protein Refseq	NP_078784
UniProt ID	Q9BTP6
Function	DNA binding; metal ion binding