



ZBED3 blocking peptide (CDBP6452)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZBED3 zinc finger, BED-type containing 3 [Homo sapiens (human)]
Official Symbol	ZBED3
Synonyms	ZBED3; zinc finger, BED-type containing 3; zinc finger BED domain-containing protein 3; 2610005H11Rik; axin-interacting protein; zinc finger, BED domain containing 3
Entrez Gene ID	84327
mRNA Refseq	NM_032367
Protein Refseq	NP_115743
UniProt ID	Q96IU2
Function	DNA binding; metal ion binding