



ZNF350 blocking peptide (CDBP6454)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZNF350 zinc finger protein 350 [Homo sapiens (human)]
Official Symbol	ZNF350
Synonyms	ZNF350; zinc finger protein 350; ZFQR; ZBRK1; zinc finger protein ZBRK1; zinc-finger protein ZBRK1; KRAB zinc finger protein ZFQR; zinc finger and BRCA1-interacting protein with a KRAB domain 1
Entrez Gene ID	59348
mRNA Refseq	NM_021632
Protein Refseq	NP_067645
UniProt ID	Q9GZX5

Pathway	Gene Expression; Generic Transcription Pathway
Function	DNA binding; metal ion binding; protein binding