



ZBTB1 blocking peptide (CDBP6455)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZBTB1 zinc finger and BTB domain containing 1 [Homo sapiens (human)]
Official Symbol	ZBTB1
Synonyms	ZBTB1; zinc finger and BTB domain containing 1; ZNF909; zinc finger and BTB domain-containing protein 1
Entrez Gene ID	22890
mRNA Refseq	NM_001123329
Protein Refseq	NP_001116801
UniProt ID	Q9Y2K1
Function	DNA binding; metal ion binding; protein heterodimerization activity; protein homodimerization

