



ZBTB3 blocking peptide (CDBP6457)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZBTB3 zinc finger and BTB domain containing 3 [Homo sapiens (human)]
Official Symbol	ZBTB3
Synonyms	ZBTB3; zinc finger and BTB domain containing 3; zinc finger and BTB domain-containing protein 3
Entrez Gene ID	79842
mRNA Refseq	NM_024784
Protein Refseq	NP_079060
UniProt ID	Q9H5J0
Function	DNA binding; metal ion binding