



ZBTB4 blocking peptide (CDBP6458)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZBTB4 zinc finger and BTB domain containing 4 [Homo sapiens (human)]
Official Symbol	ZBTB4
Synonyms	ZBTB4; zinc finger and BTB domain containing 4; ZNF903; KAISO-L1; zinc finger and BTB domain-containing protein 4; KAISO-like zinc finger protein 1
Entrez Gene ID	57659
mRNA Refseq	NM_001128833
Protein Refseq	NP_001122305
UniProt ID	Q9P1Z0
Function	metal ion binding; methyl-CpG binding; methyl-CpNpG binding; protein binding; protein

homodimerization activity; protein kinase binding; sequence-specific DNA binding
