



ZBTB5 blocking peptide (CDBP6459)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZBTB5 zinc finger and BTB domain containing 5 [Homo sapiens (human)]
Official Symbol	ZBTB5
Synonyms	ZBTB5; zinc finger and BTB domain containing 5; zinc finger and BTB domain-containing protein 5
Entrez Gene ID	9925
mRNA Refseq	NM_014872
Protein Refseq	NP_055687
UniProt ID	O15062
Function	DNA binding; metal ion binding