



ZC3H12C blocking peptide (CDBP6465)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZC3H12C zinc finger CCCH-type containing 12C [Homo sapiens (human)]
Official Symbol	ZC3H12C
Synonyms	ZC3H12C; zinc finger CCCH-type containing 12C; MCP1P3; probable ribonuclease ZC3H12C; MCP induced protein 3; MCP-induced protein 3; zinc finger CCCH domain-containing protein 12C
Entrez Gene ID	85463
mRNA Refseq	NM_033390
Protein Refseq	NP_203748
UniProt ID	Q9C0D7

Function

endonuclease activity; metal ion binding; molecular_function
