



ZC3H12D blocking peptide (CDBP6466)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZC3H12D zinc finger CCCH-type containing 12D [Homo sapiens (human)]
Official Symbol	ZC3H12D
Synonyms	ZC3H12D; zinc finger CCCH-type containing 12D; TFL; p34; MCP4; C6orf95; dJ281H8.1; probable ribonuclease ZC3H12D; tumor suppressor TFL; MCP-induced protein 4; transformed follicular lymphoma; Zinc finger CCCH domain-containing protein 12D
Entrez Gene ID	340152
mRNA Refseq	NM_207360
Protein Refseq	NP_997243
UniProt ID	A2A288

Function	endonuclease activity; metal ion binding
-----------------	--
