



## Slc39a10 blocking peptide (CDBP6476)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

### GENE INFORMATION

Gene Name	<a href="#">Slc39a10 solute carrier family 39 (zinc transporter), member 10 [ Rattus norvegicus (Norway rat) ]</a>
Official Symbol	Slc39a10
Synonyms	Slc39a10; solute carrier family 39 (zinc transporter), member 10; zinc transporter ZIP10
Entrez Gene ID	<a href="#">363229</a>
Pathway	Metal ion SLC transporters; SLC-mediated transmembrane transport; Transmembrane transport of small molecules; Transport of glucose and other sugars; Zinc influx into cells by the SLC39 gene family; Zinc transporters
Function	metal ion transmembrane transporter activity