



# SLC39A11 blocking peptide (CDBP6477)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

## GENE INFORMATION

Gene Name	<a href="#">SLC39A11 solute carrier family 39, member 11 [ Homo sapiens (human) ]</a>
Official Symbol	SLC39A11
Synonyms	SLC39A11; solute carrier family 39, member 11; ZIP11; C17orf26; zinc transporter ZIP11; ZIP-11; Zrt- and Irt-like protein 11; solute carrier family 39 (metal ion transporter), member 11
Entrez Gene ID	<a href="#">201266</a>
mRNA Refseq	<a href="#">NM_001159770</a>
Protein Refseq	<a href="#">NP_001153242</a>
UniProt ID	Q8N1S5
Function	metal ion transmembrane transporter activity