



SLC39A3 blocking peptide (CDBP6482)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SLC39A3 solute carrier family 39 (zinc transporter), member 3 [Homo sapiens (human)]
Official Symbol	SLC39A3
Synonyms	SLC39A3; solute carrier family 39 (zinc transporter), member 3; ZIP3; ZIP-3; zinc transporter ZIP3; zrt- and Irt-like protein 3; solute carrier family 39 member 3
Entrez Gene ID	29985
mRNA Refseq	NM_144564
Protein Refseq	NP_653165
UniProt ID	Q9BRY0
Pathway	Metal ion SLC transporters; SLC-mediated transmembrane transport; Senescence and

Autophagy; Transmembrane transport of small molecules; Transport of glucose and other sugars; Zinc influx into cells by the SLC39 gene family; Zinc transporters

Function	zinc ion transmembrane transporter activity
-----------------	---
