



SLC39A9 blocking peptide (CDBP6488)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SLC39A9 solute carrier family 39, member 9 [Homo sapiens (human)]
Official Symbol	SLC39A9
Synonyms	SLC39A9; solute carrier family 39, member 9; ZIP9; ZIP-9; zinc transporter ZIP9; zinc transporter SLC39A9; zrt- and Irt-like protein 9; solute carrier family 39 (zinc transporter), member 9
Entrez Gene ID	55334
mRNA Refseq	NM_001252148
Protein Refseq	NP_001239077
UniProt ID	Q9NUM3

Function

metal ion transmembrane transporter activity
