



ZMYM1 blocking peptide (CDBP6491)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZMYM1 zinc finger, MYM-type 1 [Homo sapiens (human)]
Official Symbol	ZMYM1
Synonyms	ZMYM1; zinc finger, MYM-type 1; MYM; zinc finger MYM-type protein 1; zinc finger, MYM domain containing 1
Entrez Gene ID	79830
mRNA Refseq	NM_001289088
Protein Refseq	NP_001276017
UniProt ID	Q5SVZ6
Function	nucleic acid binding; protein dimerization activity; zinc ion binding