



ZMYM4 blocking peptide (CDBP6494)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZMYM4 zinc finger, MYM-type 4 [Homo sapiens (human)]
Official Symbol	ZMYM4
Synonyms	ZMYM4; zinc finger, MYM-type 4; MYM; CDIR; ZNF262; ZNF198L3; zinc finger MYM-type protein 4; zinc finger protein 262; cell death inhibiting RNA
Entrez Gene ID	9202
mRNA Refseq	NM_005095
Protein Refseq	NP_005086
UniProt ID	Q5VZL5
Function	DNA binding; zinc ion binding