



ZNF219 blocking peptide (CDBP6495)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZNF219 zinc finger protein 219 [Homo sapiens (human)]
Official Symbol	ZNF219
Synonyms	ZNF219; zinc finger protein 219; ZFP219
Entrez Gene ID	51222
mRNA Refseq	NM_001101672
Protein Refseq	NP_001095142
UniProt ID	Q9P2Y4
Function	DNA binding; histamine receptor activity; metal ion binding; protein binding; sequence-specific DNA binding transcription factor activity