



ZNF521 blocking peptide (CDBP6497)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZNF521 zinc finger protein 521 [Homo sapiens (human)]
Official Symbol	ZNF521
Synonyms	ZNF521; zinc finger protein 521; EHZF; Evi3; LYST-interacting protein 3; early hematopoietic zinc finger protein
Entrez Gene ID	25925
mRNA Refseq	NM_015461
Protein Refseq	NP_056276
UniProt ID	Q96K83
Function	DNA binding; metal ion binding; protein domain specific binding