



ZNF536 blocking peptide (CDBP6499)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZNF536 zinc finger protein 536 [Homo sapiens (human)]
Official Symbol	ZNF536
Synonyms	ZNF536; zinc finger protein 536
Entrez Gene ID	9745
mRNA Refseq	NM_014717
Protein Refseq	NP_055532
UniProt ID	O15090
Function	metal ion binding; retinoic acid-responsive element binding