



ZNF687 blocking peptide (CDBP6500)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ZNF687 zinc finger protein 687 [Homo sapiens (human)]
Official Symbol	ZNF687
Synonyms	ZNF687; zinc finger protein 687
Entrez Gene ID	57592
mRNA Refseq	NM_020832
Protein Refseq	NP_065883
UniProt ID	Q8N1G0
Function	DNA binding; metal ion binding