



Anti-Mouse IgG2a monoclonal antibody [APC] (CDMPY-16029)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reacts with immunoglobulin G2a (IgG2a) in all tested mouse haplotype (Igh-a and b). It does not react with other isotypes.
Immunogen	Mouse Ig cocktail
Isotype	IgG
Source/Host	Rat
Species Reactivity	Mouse
Purification	Affinity chromatography
Conjugate	APC
Applications	FC
Format	Liquid
Concentration	0.2 mg/ml
Size	100 µg
Buffer	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.
Preservative	0.09% Sodium Azide
Storage	The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. Do not freeze.