



User's Manual

# Human Oxidized LDL ELISA Kit (CML-LDL)



CDN-F1761



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Human Oxidized LDL ELISA Kits (CML-LDL) are enzyme immunoassays developed for the detection and quantitation of human oxLDL in plasma, serum or other biological fluid samples.

### Storage

Upon receipt, aliquot and store the Blocking Reagent and OxLDL Standard at -20°C to avoid multiple freeze/thaw cycles. Store all other components at 4°C.

### Sensitivity

150 ng/mL