



Pseudotyped VSV-SARS-CoV-2 S-ΔG-Luciferase (Omicron variant, KP.2) (COV-PSV46)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

Pseudotyped VSV-SARS-CoV-2 S-ΔG-Luciferase (Omicron variant, KP.2) (JN.1.11.1.2) encodes the antigenomic-sense (or positive-sense) RNA of a replicaton-restricted recombinant vesicular stomatitis virus (rVSV) in which the glycoprotein (G) gene has been replaced with SARS-CoV-2 spike protein (T19I, R21T, Δ24-26, A27S, Δ69-70 V127F, G142D, Δ145, F157S, R158G, Δ211, L212I, V213G, L216F, H245N, A264D, I332V, G339H, R346T, K356T, S371F, S373P, S375F, T376A, R403K, D405N, R408S, K417N, N440K, V445H, G446S, N450D, L452W, L455S, F456L, N460K, S477N, T478K, N481K, Δ483, E484K, F486P, Q493R, Q498R, N501Y, Y505H, E554K, A570V, D614G, P621S, H655Y, 670V, N679K, P681R, N764K, D796Y, S939F, Q954H, N969K, P1143L). Because the infectivity of Pseudotyped VSV-SARS-CoV-2 S-ΔG-Luciferase is restricted to a single round of replication, the pseudotypes can be handled using BSL-2 containment practices. The pseudotype VSV particles encode Luciferase together with the VSV nucleocapsid (N), phosphoprotein (P), glycoprotein (G), and large polymerase subunit (L) in their pVSV-ΔG vector. When the VSV pseudovirus infect the target cells, Luciferase expression is proportional to the number of cells that were infected.

Specificity	SARS-CoV-2
Species	Virus
Size	1 mL
Storage	Store at -80°C. Multiple freeze/thaw cycles not recommended. When using the virus, transfer the virus from the -80° C refrigerator and melt it in an ice bath.
Ship	Frozen on dry ice

BACKGROUND

Keywords

SARS-CoV-2; Coronavirus; SARS-CoV; SARS-CoV-2 Envelope Protein; SARS-CoV-2 Envelope glycoprotein; SARS-CoV-2 glycoprotein; Severe Acute Respiratory Syndrome; Severe Acute Respiratory Syndrome-2; Severe Acute Respiratory Syndrome glycoprotein; Severe Acute Respiratory Syndrome-2 glycoprotein; SARS-CoV-2 Pseudovirus; SARS-CoV Pseudovirus; SARS Pseudovirus; JN.1; BA.2.86.1.1; BA.2.86; KP.2
