



Anti-CATHL3 polyclonal antibody (CPBT-29201RC)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Chicken CATHL3.
Specificity	Expressed in bone marrow and testis and neutrophils.
Immunogen	Synthetic peptide from Chicken Cathelicidin conjugated to an immunogenic carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Chicken
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Sequence Similarities	Belongs to the cathelicidin family.
Cellular Localization	Secreted.
Format	Liquid
Size	100 µl
Buffer	Preservative: None Constituents: Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume

of 40% for extra stability). Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	CATHL3 cathelicidin antimicrobial peptide [Gallus gallus]
Official Symbol	CATHL3
Synonyms	CATHL3; cathelicidin antimicrobial peptide; cathelicidin-1; 18 kDa cationic antimicrobial protein; Antibacterial protein FALL-39; Antibacterial protein LL-37; CAMP; CAMP_HUMAN; CAP 18; CAP-18; CAP18; Cathelicidin antimicrobial peptide; Cathelin-like protein; Cathelin-related antimicrobial peptide; CATHL3; Cationic antimicrobial protein, 18-KD; CLP; Cnlp; Cramp; CRAMP, mouse, homolog of; FALL 39; FALL-39 peptide antibiotic; FALL39; hCAP 18; hCAP-18; hCAP18; HSD26; LL37; MCLP; Peptide antibiotic, PR-39, porcine, homolog of; CATH-1; CATH-3; fowlicidin 1; fowlicidin-1; fowlicidin-3; cathelicidin-3; CAMP; CATHL1;
Entrez Gene ID	414337
Protein Refseq	NP_001001605
UniProt ID	C4PFJ7