



Anti-CHI3L3 polyclonal antibody (CPBT-30197RM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Mouse CHI3L3.
Immunogen	Synthetic peptide conjugated to KLH derived from within residues 350 to the C-terminus of Human Chitinase 3 like protein 3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Secreted
Format	Liquid
Size	100 µg
Buffer	Preservative: 0.02% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	Chi3l3 chitinase 3-like 3 [Mus musculus]
Official Symbol	CHI3L3
Synonyms	CHI3L3; chitinase 3-like 3; chitinase-3-like protein 3; secreted protein Ym1; eosinophil chemotactic cytokine; eosinophil chemotactic factor-L; Ym1; ECF-L; AI505981;
Entrez Gene ID	12655
Protein Refseq	NP_034022
UniProt ID	O35744
Function	beta-N-acetylhexosaminidase activity; catalytic activity; cation binding; chitin binding; NOT chitinase activity; hydrolase activity, hydrolyzing O-glycosyl compounds; sugar binding;