



# Rabbit Anti-Human COL20A1 Polyclonal Antibody (CPBT-31628RH)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant Protein, antigen sequence: LVSAIYAAGRSEAVSATGQTACPALRPDGS LPGFDLMVAFSLVEKAYASIRGVAMEPSAF GGTPTFTLFKDAQLTRRVSDVYPA
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Antigen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	40% glycerol and PBS (pH 7.2).
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at +4°C for short term storage. Long time storage is recommended at -20°C. Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

## BACKGROUND

<b>Introduction</b>	COL20A1 is a probable collagen protein.
<b>Keywords</b>	COL20A1; collagen, type XX, alpha 1; bA261N11.4; collagen alpha-1(XX) chain; collagen-like protein;

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">57642</a>
<b>UniProt ID</b>	<a href="#">Q9P218</a>