



Rabbit anti-Saccharomyces EAF3 Polyclonal antibody (CPBT-32848RS)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Saccharomyces EAF3.
Immunogen	Recombinant full length protein of <i>S. cerevisiae</i> Eaf3.
Isotype	IgG
Source/Host	Rabbit
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ChIP, IP, WB
Cellular Localization	Nuclear
Format	Liquid
Size	100 µl
Buffer	Preservative: 0.02% Sodium Azide
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [EAF3 Eaf3p \[Saccharomyces cerevisiae S288c \]](#)

Official Symbol	EAF3
Synonyms	EAF3; Eaf3p; Chromatin modification related protein EAF3; Eaf 3; Eaf3p; ESA1 associated factor 3; Esa1p-associated factor, nonessential component of the NuA4 acetyltransferase complex, homologous to Drosophila dosage compensation protein MSL3; plays a role in regulating Ty1 transposition;
Entrez Gene ID	856134
Protein Refseq	NP_015348
UniProt ID	Q12432
Function	histone acetyltransferase activity;