



Anti-FNIP1 (N-terminal) polyclonal antibody (CPBT-33934RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Human FNIP1.
Immunogen	A peptide at the N-terminus of FNIP1, MAPTLFQKLFSKRTGLGAC-amide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	50 µl
Buffer	Preservative: 0.1% Sodium Azide
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [FNIP1 folliculin interacting protein 1 \[Homo sapiens \]](#)

Official Symbol	FNIP1
Synonyms	FNIP1; folliculin interacting protein 1; folliculin-interacting protein 1; KIAA1961; FLCN-interacting protein; folliculin interacting protein 1;
Entrez Gene ID	96459
Protein Refseq	NP_001008738
UniProt ID	Q8TF40
Chromosome Location	5q23.3
Function	protein binding;