



Anti-FIGNL1 (aa 1-50) polyclonal antibody (CPBT-34631RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Human FIGNL1.
Immunogen	Synthetic peptide, corresponding to a region between amino acids 1-50 of Human FIGNL1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µl
Buffer	Preservative: 0.09% Sodium Azide Constituents: 8mM PBS, 60mM Citrate, 150mM Tris, pH 7 to 8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

GENE INFORMATION

Gene Name [FIGNL1 fidgetin-like 1 \[Homo sapiens \]](#)

Official Symbol	FIGNL1
Synonyms	FIGNL1; fidgetin-like 1; fidgetin-like protein 1; Fidgetin like protein 1;
Entrez Gene ID	63979
Protein Refseq	NP_071399
UniProt ID	Q6PIW4
Chromosome Location	7p12.2
Function	ATP binding; hydrolase activity; magnesium ion binding; nucleoside-triphosphatase activity; nucleotide binding;