



Anti-IFIT5 (aa 79-294) polyclonal antibody (CPBT-37313RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Human IFIT5.
Immunogen	Recombinant fragment, corresponding to amino acids 79-294 of Human IFIT5 (BC025786).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Sequence Similarities	Belongs to the IFIT family. Contains 8 TPR repeats.
Format	Liquid
Size	100 µg
Buffer	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 0.5% BSA, 99% PBS
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	IFIT5 interferon-induced protein with tetratricopeptide repeats 5 [Homo sapiens]
Official Symbol	IFIT5
Synonyms	IFIT5; interferon-induced protein with tetratricopeptide repeats 5; retinoic acid and interferon inducible protein (58kD); RI58; FLJ53857; FLJ92678; IFIT-5; IFIT5; IFIT5_HUMAN; Interferon induced protein with tetratricopeptide repeats 5; Interferon-induced protein with tetratricopeptide repeats 5; Retinoic acid and interferon inducible 58 kDa protein; Retinoic acid- and interferon-inducible 58 kDa protein; RI58; IFIT-5; retinoic acid- and interferon-inducible 58 kDa protein; retinoic acid- and interferon-inducible protein (58kD);
Entrez Gene ID	24138
Protein Refseq	NP_036552
UniProt ID	Q13325
Chromosome Location	10q23.31
Function	binding; molecular_function;