



Anti-LAP3 (full length) polyclonal antibody (CPBT-37379MH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Polyclonal antibody to Human LAP3.
Immunogen	LAP3full-length human protein.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic
Format	Liquid
Size	50 µg
Buffer	Preservative: NoneConstituents: 1X PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	LAP3 leucine aminopeptidase 3 [Homo sapiens]
Official Symbol	LAP3
Synonyms	LAP3; leucine aminopeptidase 3; PEPS,peptidase S; cytosol aminopeptidase; LAP; LAPEP; Cytosol aminopeptidase; LAP 3; LAP; LAPEP; Leucine aminopeptidase 3; Leucyl aminopeptidase; PEPS; Peptidase S; Proline aminopeptidase; Prolyl aminopeptidase; LAP-3; peptidase S; leucyl aminopeptidase; prolyl aminopeptidase; proline aminopeptidase; PEPS;
Entrez Gene ID	51056
Protein Refseq	NP_056991
UniProt ID	P28838
Chromosome Location	4p15.33
Pathway	Arginine and proline metabolism, organism-specific biosystem; Arginine and proline metabolism, conserved biosystem; Glutathione metabolism, organism-specific biosystem; Glutathione metabolism, conserved biosystem; Metabolic pathways, organism-specific biosystem;
Function	aminopeptidase activity; magnesium ion binding; manganese ion binding; metalloexopeptidase activity; peptidase activity; zinc ion binding;