



Anti-LRMP (C-terminal) polyclonal antibody (CPBT-38057RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Human LRMP.
Antigen Description	The protein encode dby this gene is expressed in a developmentally regulated manner in lymphoid cell lines and tissues. The protein is localized to the cytoplasmic face of the endoplasmic reticulum.
Specificity	Expressed at high levels in pre B-cells, mature B-cells and pre T-cells. Expressed at low levels in mature T-cells and plasma B-cells.
Immunogen	Synthetic peptide at the C terminal of human JAW1 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Cellular Localization	Cytoplasm and Endoplasmic reticulum membrane.
Format	Liquid
Size	200 µl
Buffer	Preservative: 0.02% Sodium AzideConstituents: 2% BSA

Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	LRMP lymphoid-restricted membrane protein [Homo sapiens]
Official Symbol	LRMP
Synonyms	LRMP; lymphoid-restricted membrane protein; JAW1; Lrmp; LRMP_HUMAN; Lymphoid restricted membrane protein; Processed lymphoid-restricted membrane protein; Protein Jaw1; protein Jaw1;
Entrez Gene ID	4033
Protein Refseq	NP_001191055
UniProt ID	Q12912
Chromosome Location	12p12.1