



Anti-LYPLAL1 (aa 87-127) polyclonal antibody (CPBT-39198RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Human LYPLAL1.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 87-127 of Human LYPLAL1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Sequence Similarities	Belongs to the AB hydrolase 2 family.
Cellular Localization	Cytoplasm.
Format	Liquid
Size	100 µg
Buffer	Preservative:0.09% Sodium azideConstituent:99% PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name	LYPLAL1 lysophospholipase-like 1 [Homo sapiens]
Official Symbol	LYPLAL1
Synonyms	LYPLAL1; lysophospholipase-like 1; lysophospholipase-like protein 1; Q96AV0; BC027340; FLJ99730; KIAA1238; LYPL1_HUMAN; LYPLAL1; Lysophospholipase like 1; Lysophospholipase-like protein 1; MGC28394; OTTHUMP00000035229; OTTHUMP00000035229; OTTHUMP00000035230;
Entrez Gene ID	127018
Protein Refseq	NP_620149
UniProt ID	Q5VWZ2
Chromosome Location	1q41
Function	hydrolase activity; lysophospholipase activity;