



# Anti-PAQR7 (internal region) polyclonal antibody (CPBT-39421RH)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Polyclonal antibody to Human PAQR7.
<b>Immunogen</b>	A synthesized peptide derived from an internal region of human MPRA.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, WB
<b>Cellular Localization</b>	Membrane; Multi pass membrane protein
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), 150mM Sodium chloride, pH 7.4
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">PAQR7 progestin and adipoQ receptor family member VII [ Homo sapiens ]</a>
<b>Official Symbol</b>	PAQR7
<b>Synonyms</b>	PAQR7; progestin and adipoQ receptor family member VII; membrane progestin receptor alpha; MPRA; mSR; Membrane progestin receptor alpha; mPR alpha; MRPA; mSR; PAQR7; PGLP; Progestin and adipoQ receptor family member VII; mPR alpha; 2310021M12Rik; progestin and adipoQ receptor family member 7; PGLP;
<b>Entrez Gene ID</b>	<a href="#">164091</a>
<b>Protein Refseq</b>	<a href="#">NP_848509</a>
<b>UniProt ID</b>	<a href="#">Q86WK9</a>
<b>Chromosome Location</b>	1p35.3
<b>Function</b>	receptor activity; steroid binding;