



Anti-FAM132A (aa 40-139) polyclonal antibody (CPBT-42120MM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Polyclonal antibody to Mouse FAM132A.
Immunogen	Fusion protein: ITTSNREGFPGSVKPPEASGPELSDAHMTWLNFVRRPDDG SSRKRCRGRDKKSRGLSGLPGPPGPPGPPGPPGSPGVGT PEALLQEFQEILKEATELRF, corresponding to amino acids 40/139 of Mouse RAI150
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Mouse
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cellular component: extracellular space
Format	Liquid
Size	100 µl
Buffer	Constituents: 50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	Fam132a family with sequence similarity 132, member A [Mus musculus]
Official Symbol	FAM132A
Synonyms	FAM132A; family with sequence similarity 132, member A; protein FAM132A; C1qTNF12; hypothetical protein LOC67389; RIKEN cDNA 1110035L05; adipolin; C1q domain containing 2; C1qdc2; 1110035L05Rik;
Entrez Gene ID	67389
Protein Refseq	NP_080401
UniProt ID	Q8R2Z0
Function	molecular_function;