



Anti-PfRh5 (aa 45-144) polyclonal antibody (CPBT-43246MP)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Polyclonal antibody to Plasmodium PfRh5.
Immunogen	Fusion protein: KSTEEKDDIKNGKDIKKEIDNDKENIKTNNAKDHSTYIK SYLNTNVNDGLKYLFIPIHNSFIKKYSVFNQINDGMLLNE KNDVKNNEDYKNVDYKKNVNF, corresponding to amino acids 45/144 of Plasmodium falciparum PFD1145c
Isotype	IgG
Source/Host	Mouse
Species Reactivity	P. falciparum
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Constituents: 50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	PfRh5 reticulocyte-binding protein homologue 5 [Plasmodium falciparum 3D7]
Official Symbol	PfRh5
Synonyms	PfRh5; reticulocyte-binding protein homologue 5; Hypothetical protein PFD1145c;
Entrez Gene ID	812437
UniProt ID	Q8IFM5