



Rabbit Anti-Human RAB32 Polyclonal Antibody (CPBT-44065MH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: NKDSSQSPSQVDQFCKEHGFAGWFETSAKDNINIEEAARFLVEKILVNHQSFPNEENDVD KIKLDQETLRAENKSQCC
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purified
Conjugate	Unconjugated
Applications	IHC, WB
Format	Liquid
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2).
Preservative	0.02% Sodium Azide
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C. Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

BACKGROUND

Introduction	Small GTP-binding proteins of the RAB family, such as RAB32, play essential roles in vesicle and granule targeting (Bao et al., 2002 [PubMed 11784320]).[supplied by OMIM, Aug 2009]
Keywords	RAB32; RAB32, member RAS oncogene family; ras-related protein Rab-32;

GENE INFORMATION

Entrez Gene ID	10981
UniProt ID	Q13637