



Rabbit Anti-Human SMC1B Polyclonal Antibody (CPBT-45544RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1012-1161 of human SMC1B (NP_001278430.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB
Positive Control	Rat testis
Format	Liquid
Size	50 µl, 100 µl
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords SMC1B; structural maintenance of chromosomes 1B; SMC1L2; SMC1BETA; structural

maintenance of chromosomes protein 1B; SMC protein 1B; SMC1 structural maintenance of chromosomes 1-like 2; SMC1 (structural maintenance of chromosomes 1, yeast)-like 1; mitosis-specific chromosome segregation protein like protein beta;

GENE INFORMATION

Entrez Gene ID	27127
----------------	-----------------------

UniProt ID	Q8NDV3
------------	------------------------
