



Anti-SERAC1 (aa 555-585) polyclonal antibody (CPBT-45959RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Human SERAC1.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 555-585 of Human SERAC1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Sequence Similarities	Belongs to the SERAC1 family.
Cellular Localization	Membrane.
Format	Liquid
Size	100 µg
Buffer	Preservative: 0.09% Sodium AzideConstituents: PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name	SERAC1 serine active site containing 1 [Homo sapiens]
Official Symbol	SERAC1
Synonyms	SERAC1; serine active site containing 1; protein SERAC1; FLJ14917; FLJ30544; 4930511N22Rik; D17Ert141e; FLJ14917; Protein SERAC1; RGD1565064; RP1-257C19.1; Serac1; Serine active site containing 1; SRAC1_HUMAN;
Entrez Gene ID	84947
Protein Refseq	NP_116250
UniProt ID	Q96JX3
Chromosome Location	6q25.3
Function	binding; hydrolase activity, acting on ester bonds;