



Anti-SLTM (full length) polyclonal antibody (CPBT-46557MH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Polyclonal antibody to Human SLTM.
Immunogen	Full length human SLTM protein (ENSP00000369883, 1 a.a. ~ 478 a.a)
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein G purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear. Note: Detected in punctate structures.
Format	Liquid
Size	50 µg
Buffer	Preservative: None Constituents: 1X PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	SLTM SAFB-like, transcription modulator [Homo sapiens]
Official Symbol	SLTM
Synonyms	SLTM; SAFB-like, transcription modulator; SAFB-like transcription modulator; FLJ13213; Met; DKFZp762G052; FLJ10005; FLJ13213; MET; Modulator of estrogen induced transcription; SAFB like transcription modulator; modulator of estrogen-induced transcription;
Entrez Gene ID	79811
Protein Refseq	NP_001013865
UniProt ID	Q9NWH9
Chromosome Location	15q21.3
Function	RNA binding; nucleotide binding;