



## Anti-STRA8 (aa 1-100) polyclonal antibody (CPBT-47319RM)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Polyclonal antibody to Mouse STRA8.
<b>Immunogen</b>	Synthetic peptide derived from within residues 1 - 100 of Mouse Stra8. ( Immunogen available as <a href="#">DAG-P1203</a> )
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse, Rat, Hamster
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P, WB, IHC-Fr, ICC/IF, IHC-FoFr
<b>Cellular Localization</b>	Cytoplasmic
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	Preservative: 0.02% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

### GENE INFORMATION

<b>Gene Name</b>	<a href="#">Stra8 stimulated by retinoic acid gene 8 [ Mus musculus ]</a>
<b>Official Symbol</b>	STRA8
<b>Synonyms</b>	STRA8; stimulated; stimulated by retinoic acid gene 8 protein; Stimulated by retinoic acid Gene 8;
<b>Entrez Gene ID</b>	<a href="#">20899</a>
<b>Protein Refseq</b>	<a href="#">NP_033318</a>
<b>UniProt ID</b>	<a href="#">P70278</a>