



Rabbit Anti-Human KIAA0922 Polyclonal Antibody (CPBT-48021RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: LNCTLENGVPCVIQESAPVHNSFIDWSATCEGQFSSAYCPLELNDYNAPPEENMNYANGF PCPADVQTDFFIDHNSQSTWNTPPNMPAAWGHASFISSP
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purified
Conjugate	Unconjugated
Applications	IHC, ICC-IF
Format	Liquid
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2).
Preservative	0.02% Sodium Azide
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C. Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

BACKGROUND

Introduction	TMEM131L belongs to the TMEM131 family. There are 4 isoforms. The exact function of this protein is not known.
---------------------	--

Keywords	KIAA0922; TMEM131L; transmembrane protein 131-like;
-----------------	---

GENE INFORMATION

Entrez Gene ID	23240
-----------------------	-----------------------

UniProt ID	A2VDJ0
-------------------	------------------------
