



Anti-KIAA0922 (aa 1275-1305) polyclonal antibody (CPBT-48022RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Human KIAA0922.
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 1275-1305 of Human TMEM131L, conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane; Single pass type I membrane protein
Format	Liquid
Size	100 µg
Buffer	Preservative: 0.09% Sodium AzideConstituents: PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name	KIAA0922 KIAA0922 [Homo sapiens]
Official Symbol	KIAA0922
Synonyms	KIAA0922; transmembrane protein 131-like; DKFZp586H1322; TMEM131L; KIAA0922; Transmembrane protein 131 like;
Entrez Gene ID	23240
Protein Refseq	NP_001124479
UniProt ID	A2VDJ0
Chromosome Location	4q31.3