



Rabbit Anti-Human UGT3A2 Polyclonal antibody (CPBT-48305RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	UGT3A2 fusion protein, sequence: QVISWLAPEDHQREFKKSFDFFLEETLGGRGKFENLLNVLEYLALQCSHFLNRKD (75-129 aa encoded by BC103925)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, ELISA
Positive Control	human testis tissue
Format	Liquid
Size	50 µl, 100 µl
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquoting is unnecessary for -20°C storage.

BACKGROUND

Keywords UGT3A2; UDP glycosyltransferase 3 family, polypeptide A2; UDP-glucuronosyltransferase 3A2;
UDPGT 3A2;

GENE INFORMATION

Entrez Gene ID [167127](#)

UniProt ID [Q3SY77](#)
