



Anti-UGT3A2 (aa 384-413) polyclonal antibody (CPBT-48306RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Human UGT3A2.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 384-413 of Human UGT3A2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane; Single-pass type I membrane protein
Format	Liquid
Size	100 µg
Buffer	Preservative: 0.09% Sodium AzideConstituents: PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name	UGT3A2 UDP glycosyltransferase 3 family, polypeptide A2 [Homo sapiens]
Official Symbol	UGT3A2
Synonyms	UGT3A2; UDP glycosyltransferase 3 family, polypeptide A2; UDP-glucuronosyltransferase 3A2; EC 2.4.1.17; MGC119426; MGC119429; OTTHUMP00000221219; UDP glucuronosyltransferase 3A2; UDP glycosyltransferase 3 family polypeptide A2; UDPGT 3A2; UDPGT 3A2; OTTHUMP00000161736; OTTHUMP00000221219;
Entrez Gene ID	167127
Protein Refseq	NP_001161788
UniProt ID	Q3SY77
Chromosome Location	5p13.2
Function	glucuronosyltransferase activity;