



Anti-TRAPPC6A (aa 42-72) polyclonal antibody (CPBT-48411RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Human TRAPPC6A.
Antigen Description	TRAPPC6A is a subunit of the TRAPP I and II complexes (Gwynn et al., 2006
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within N terminal amino acids 42-72 of Human TRAPPC6A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Cellular Localization	Golgi apparatus; cis-Golgi network. Endoplasmic reticulum.
Format	Liquid
Size	100 µg
Buffer	Preservative: 0.09% Sodium AzideConstituents: PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name	TRAPPC6A trafficking protein particle complex 6A [Homo sapiens]
Official Symbol	TRAPPC6A
Synonyms	TRAPPC6A; trafficking protein particle complex 6A; trafficking protein particle complex subunit 6A; HSPC289; MGC2650; TRS33; 4930519D19Rik; AI480686; MGC117907; mhyp; OTTMUSP00000025166; RGD1309160; TRS33; 1810073E21Rik; HSPC289; OTTMUSP00000025165; rCG_54352; Trafficking protein particle complex subunit 6A; TRAPP complex subunit 6A; TRAPPC6Adelta29-42; TRAPP complex subunit 6A;
Entrez Gene ID	79090
Protein Refseq	NP_077013
UniProt ID	B7TZ90
Chromosome Location	19q13.32
Function	binding;