



Anti-USP35 (aa 951-1000) polyclonal antibody (CPBT-49090RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Human USP35.
Specificity	Expressed in testis, pancreas and skeletal muscle.
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 951-1000 (MEAISKDNIL YLQEQEKEAR SRAAYISALP TSPHWGRGFD EDKDEDEGSP) of Human USP35.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Sequence Similarities	Belongs to the peptidase C19 family.
Format	Liquid
Size	50 µg
Buffer	Preservative: None Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze /

thaw cycles.

GENE INFORMATION

Gene Name	USP35 ubiquitin specific peptidase 35 [Homo sapiens]
Official Symbol	USP35
Synonyms	USP35; ubiquitin specific peptidase 35; ubiquitin specific protease 35; ubiquitin carboxyl-terminal hydrolase 35; KIAA1372; USP34; Deubiquitinating enzyme 35; Gm1088; Gm493; KIAA1372; Ubiquitin carboxyl terminal hydrolase 35; Ubiquitin carboxyl-terminal hydrolase 35; Ubiquitin specific peptidase 35; Ubiquitin specific processing protease 35; Ubiquitin specific protease 35; Ubiquitin thioesterase 35; Ubiquitin thiolesterase 35; Ubiquitin-specific-processing protease 35; UBP35_HUMAN; USP35; ubiquitin thioesterase 35; deubiquitinating enzyme 35; ubiquitin thiolesterase 35; ubiquitin-specific-processing protease 35;
Entrez Gene ID	57558
Protein Refseq	NP_065849
UniProt ID	B3KPE1
Chromosome Location	11q13.4
Function	binding; cysteine-type peptidase activity; peptidase activity; ubiquitin thiolesterase activity;