



Anti-KIAA1033 (internal region) polyclonal antibody (CPBT-49784RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Human KIAA1033.
Immunogen	Synthetic peptide, derived from an internal domain of Human WASH complex subunit 7.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P
Sequence Similarities	Belongs to the WASHS7 family.
Format	Liquid
Size	200 µl
Buffer	Constituent:Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	KIAA1033 KIAA1033 [Homo sapiens]
Official Symbol	KIAA1033
Synonyms	KIAA1033; WASH complex subunit 7; Kiaa1033; WAHS7_HUMAN; WASH complex subunit 7;
Entrez Gene ID	23325
Protein Refseq	NP_056090
UniProt ID	A0A087X256
Chromosome Location	12q24.11