



Anti-ZDHHC21 (aa 68-98) polyclonal antibody (CPBT-50179RH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Human ZDHHC21.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within N terminal amino acids 68-98 of Human ZDHHC21.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane; Multi-pass membrane protein
Format	Liquid
Size	100 µg
Buffer	Preservative: 0.09% Sodium AzideConstituents: PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name	ZDHHC21 zinc finger, DHHC-type containing 21 [Homo sapiens]
Official Symbol	ZDHHC21
Synonyms	ZDHHC21; zinc finger, DHHC-type containing 21; probable palmitoyltransferase ZDHHC21; DNZ1; HSPC097; 9130404H11Rik; AL024349; D130004H04Rik; DHHC 21; DNZ1; GABA A receptor-associated membrane protein 3; HSPC097; OTTMUSP00000000336; OTTMUSP00000022012; Probable palmitoyltransferase ZDHHC21; RP23-76J17.1; Zinc finger DHHC domain containing protein 21; zinc finger, DHHC domain containing 21; zinc finger DHHC domain-containing protein 21; DHHC-21; 9130404H11Rik;
Entrez Gene ID	340481
Protein Refseq	NP_848661
UniProt ID	Q8IVQ6
Chromosome Location	9p22.3
Pathway	Metabolism, organism-specific biosystem; Metabolism of nitric oxide, organism-specific biosystem; eNOS activation, organism-specific biosystem; eNOS activation and regulation, organism-specific biosystem;
Function	metal ion binding; transferase activity, transferring acyl groups; zinc ion binding;