

Anti-Rabbit IgG polyclonal antibody [Biotin] (CPBT-67182SR)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep anti Rabbit IgG antibody recognizes rabbit IgG (Heavy and light chains). Cross reactivity with bovine, human, mouse or rat serum components has been minimised by adsorption. Sheep anti Rabbit IgG antibody has been used extensivley as a secondary reagent for the detection of rabbit polyclonal antibodies in immunohistological and ELISA applications.
Specificity	IGGC
Immunogen	Rabbit IgG from pooled sera
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Rabbit
Conjugate	Biotin
Applications	IHC-Fr; ELISA; IHC-P; WB
Format	Purified IgG conjugated to Biotin - liquid
Size	1 ml
Preservative	See individual product datasheet
Storage	in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

45-1 Ramsey Road, Shirley, NY 11967, USA

Email: info@creative-diagnostics.com

Tel: 1-631-624-4882 Fax: 1-631-938-8221

Gene Name	LOC100009097 lg gamma chain C region [Oryctolagus cuniculus (rabbit)]
Official Symbol	IGGC
Synonyms	LOC100009097; IGGC; Ig gamma chain C region; Ig gamma H-chain; immunoglobulin gamma1 constant region;
Entrez Gene ID	100009097
UniProt ID	P01870
Pathway	African trypanosomiasis; Allograft rejection; Amoebiasis; Asthma;